

Denpro

ORTHODONTIC INSTRUMENTS

- Orthodontic Pliers
- Ortho Instruments
- Impression Trays
- Hollow Wares




Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



709-003-130
 Young
 13 cm, 5"

Dahtbiegezange
 Wire Bending Plier
 Pince à courber les fils
 Alicates para doblar alambre



709-006-165
 Tryfus
 16.5 cm, 6 1/2"

Kronenformzange
 Crown Forming Plier
 Pince pour former les couronnes
 Alicates para formar coronas




709-009-120
 Johnson
 12 cm, 4 1/4"


Konturzangen
 Contouring Pliers
 Pince à bouteroller
 Alicates para contornear




709-012-130
 Reynolds
 13 cm, 5"




709-015-135
 Gordon
 13.5 c, 5 1/4"



709-018-135
 USA Modell
 13.5 cm, 5 1/4"



709-021-135
 Goslee
 13.5 cm, 5 1/4"



709-024-135
 Frevert
 13 cm, 5"

Konturzangen
 Contouring Pliers
 Pinces à bouteroller
 Alicates para contornear

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



709-027-150
 Kunz
 15 cm, 6"

kronenformzange
 Crown Forming Plier
 Pince pour former les couronnes
 alicates para formar coronas



709-030-150
 Böhm
 15 cm, 6"

Teleskop-kronenzange
 Telescope Crown Plier
 Pince pour couronnes télescope
 alicates para coronas telescopio



709-033-130
 Waldron
 13 cm, 5"

Drahtbiegezange
 wire bending Plier
 Pince à courber les fils
 alicates para doblar alambre




709-036-140
 Triplex
 14 cm, 5 1/2"

konturzangen
 Contouring Pliers
 Pince à bouteroller
 alicates para contornear



709-039-150
 universal
 15 cm, 6"



709-042-150
 16 cm, 6 1/2"

Quecksilberauspreßzange
 Amalgam Squeezing Plier
 Pince à comprimer l'amalgame
 pinza para comprimir amalgama

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



709-068-130
 13 cm, 5"



709-071-130
 Schwarz
 13 cm, 5"




709-074-130
 Schwarz
 13 cm, 5"




709-077-130
 Schwarz
 14 cm, 5 1/4"


Draht- und Bügelbiegezangen
 Arch and Spring Bending Pliers
 Pinces à plier les fils et les barres
 Alicates para doblar alambre y barras




709-080-125
 12.5 cm, 4 3/4"



709-083-125
 Angle
 12.5 cm, 4 3/4"



709-086-125
 Adams
 12.5 cm, 4 3/4"




709-089-130
 13 cm, 5"

Drahtbiegezangen
 Wire Bending Pliers
 Pinces à plier les fils
 Alicates para doblar alambre


Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis




709-092-140
 14 cm, 5 1/2"

Flachzange zum Einsetzen
 von Einschraubspitzen
 Flat nose plier for screw points
 Pince plates pour pointes à vis
 Alicates de picos chatos
 para puntas a rosca




709-095-160
 16 cm, 6 1/4"




709-098-160
 16 cm, 6 1/4"

Flachzangen
 Flat nose pliers
 Pinces plates
 Alicates de picos chatos



709-101-140
 14 cm, 5 1/2"



709-104-140
 14 cm, 5 1/2"



709-107-135
 13.5 cm, 5 1/4"



709-110-135
 13.5 cm, 5 1/4"

Flachzangen
 Flat nose pliers
 Pinces plates
 Alicates de picos chatos

PLIERS FOR ORTHODONTICS & PROSTHETICS


Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis




Slim profile, tapered tips offer greater access versatility.

709-113-130
Slim Weingart Pliers
 13 cm, 5"




709-116-130
Straight How Plier
 13 cm, 5"
 All the features that makes the standard how plier (below) popular, are here, but with smaller (3/32") diameter tip.




709-119-130
Long Handle Weingart Plier
 13 cm, 5"
 This model features longer handles and a head size for a remarkably secure grip when bending wires and working in hard to reach areas.



709-122-130
Long Handle Weingart Plier
 13 cm, 5"
 A long-handled version of the most used & useful orthodontic plier. Added length means better access, greater visibility and less fatigue. Ideal for inserting wires through buccal tubes. Serrated tips securely hold wire during use.



709-125-130
Weingart Utility Plier With Offset Head
 13 cm, 5"
 45 angled head gives you easier access in difficult to reach posterior and lingual areas. Tips are serrated for secure holding. Bowed jaws eliminate lip pinching.




709-128-130
Straight How Plier
 13 cm, 5"
 one of the most widely used orthodontic hand instruments. Traditional size, bowed jaws that prevent lip pinches, with standard 1/8" tip size.

PLIERS FOR ORTHODONTICS & PROSTHETICS

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



709-131-130
Lingual Arch
Forming Pliers
 13 cm, 5"




Double back & triple back bends in .030" & .036" (.76mm & .91mm) lingual arch wires.




709-134-130
Weingart Plier
 13 cm, 5"
 This multi-purpose instrument features an offset beak for a convenient working angle, fully serrated tips, and bowed jaws that won't pinch lips.




709-137-130
How Plier
 13 cm, 5"
 Both standard and small sizes of this plier have sure-hold serrated gripping tips, and non-contacting bowed jaws that won't pinch the patient's lip. Small size for younger children.




709-140-130
Reynolds Contouring
 13 cm, 5"



709-143-130
Wire Bending With Groved
 13 cm, 5"



709-145-130
Peeso's Collar
 13 cm, 5"




709-148-130
Baade, Band & Shell Removing
 13 cm, 5"

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik


Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis




709-151-130
Tweed Loop
Forming Pliers
 13 cm, 5"




709-154-130
Three Jaw Intraoral Plier
 13 cm, 5"
 Produces superior clasps and intororal loops time after time. One beak has a center slot that accommodates the narrow center shaft on the opposing jaw. Designed for light round and rectangular wires of .030" diameter.




709-157-130
Three Jaw Lab Plier
 13 cm, 5"
 Heavy duty construction to fabricate clasps. Hawley retainers or make adjustments to lip bumpers. Longer handles and slightly larger head size provide force to quickly bend .045" diameter round and rectangular wires.



709-160-130
Adams Plier
 13 cm, 5"
 Built to take the pressure required to work .030" Adams clasps and to reduce discomfort for the user. Both beaks are flat surfaces, with jaws parallel at .040" opening. Rugged construction ensures long working life.




709-163-130
Crown/Band crimping Plier
 13 cm, 5"
 Slim tip makes it simple to crimp and reform crowns and bands for a snug fit. Slim design of the matching concave and convex beaks assures better control.



709-166-130
Lab Wire Bending Plier
 13 cm, 5"
 Heavy duty bird beak plier has a serrated rectangular tip and a smooth round tip, with jaws parallel at .040" opening. Tips are sandblasted for non-slip gripping. Slightly longer than standard size, and appreciable heavier to allow the application of greater force needed to bend .056" wires.

*Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik*

*Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis*



709-169-130
Universal Cut & Hold
Distal End Cutter
13 cm, 5"



Can easily cut all types
of wires ranging from
.012" / .30mm to .021"
x .025" (.53mm x .64mm)



709-171-130
Ligature Forming Plier
13 cm, 5"
Create your own pre-formed
ligatures with a single squeeze.
Fabricates 5mm loops with
consistent accuracy.



709-174-130
**Intraoral "V" Notch
Stop Plier**
13 cm, 5"
Now make precise 1mm
"V" stops right in the mouth.
Prevents wires from
disengaging from the buccal
tubes. Notches up to .018"
round and .016" x .022"
rectangular.



709-177-130
Loop Tie Back Plier
13 cm, 5"
Form 2mm, 3mm, 4mm
and 5mm tic back and
retraction loops in square
and rectangular wires with
two quick squeezes of the
handles.




709-180-130
Extraoral "V" Notch Stop Plier
13 cm, 5"
A single closure forms a precise 1mm "V" stop
that prevents wire from disengaging from the
buccal tube. Eliminates unscheduled
appointments to reseat wires and saves you
valuable chair time.



709-183-130
Intraoral Detailing Plier
13 cm, 5"
Fabricate delicate offsets in archwire –
accurately, quickly. Ideal for bends to engage
lingually positioned teeth. Also good for by-pass
and 2nd order bends. Makes 1mm bends.

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



709-186-130
Three Jaw Pliers
 13 cm, 5"



Useful for contouring and bending wire up to .030" / .76mm



709-189-130
Double-sided Detailing Plier
 13 cm, 5"
 Intraoral detailing without disengaging wire means significant time savings. Step up, down, in or out. Tips are at right angle from handle to allow easy access to the arch wire.

709-192-130
Optical Plier For Loops
 13 cm, 5"
 This durable plier has fine tips for making precision bends and adjustments, and gives you easy access to hard-to-reach intraoral areas. The upper jaw is rounded, fitting into a concave groove in the lower. Use wires up to .028".



709-195-130
Band Seater
 13 cm, 5"
 The serrated tips of this precision instrument fit between the tie-wings of the smallest twin brackets, giving you the secure grip needed to place bands firmly and quickly. Plier jaws contact only at tip, eliminating pinching lips.


709-198-130
Hollow Chop Plier
 13 cm, 5"
 designed to enable you to create arches in any round or rectangular stainless steel wire up to .036". The smooth rounded beaks facilitate easy contouring.

709-201-130
Contouring Ball Plier
 13 cm, 5"
 Matched ball and socket tips for accurate, convenient forming and contouring with large tips for band contouring with small tips for anterior bands and small crowns.

709-204-130
Headgear Plier - Type 1
 13 cm, 5"
 Heavy duty plier used to bend the inner wire of a facebow. Won't weaken wire by nicking or scraping. Jaw design features symmetrical faces with three grooves in each.

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



709-204-130
Bird Beak Pliers
 13 cm, 5"



Versatile loop forming pliers for round wires up to .030" / .76mm. Cone tip and opposing tip are .025"/64mm



709-207-130
Headgear Plier - Type 2
 13 cm, 5"
 This heavyweight plier safely bends the toughest outer facebows. Solid stainless steel forging, short beaks and triple head construction provide all the force you need.



709-210-130
Three Jaw Bending Plier
 13 cm, 5"
 Ideal for adjusting and contouring wires, clasps and retainers smaller than .020 at tip, and no larger than .028 at base of jaws. Bends up to .030".



709-213-130
Bird Beak Plier Without Cutter
 13 cm, 5"
 facilitates bending round or rectangular wires up to .30" diameter. This versatile instrument features hardened inserts on all bending surfaces, a sure non-slip grip, and long-lasting jaw alignment for unmatched precision.



709-216-130
Light Wire Bird Beak Plier
 13 cm, 5"
 for fabricating omega, helix or retraction loops in wires up to .020". Choose the cutter model for bending and cutting wires up to .020"



709-219-130
Secure Hold Bird Beak
 13 cm, 5"
 whether you're working with round or rectangular wire, these pliers give you a remarkably stable platform for difficult-to-fabricate clasps - - thanks to the serrated, non-slip surface of the flat beak and a closing groove.




709-222-130
Lingual Arch Forming Plier
 13 cm, 5"

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis




709-225-130
Angled Direct Bond
Bracket Remover T.c
 13 cm, 5"




709-228-130
Debonding Plier For
Ceramic & Plastic Bracket
 13 cm, 5"
 unique design features flat tips separated by 2mm gap at closure to grip securely without crushing of fracturing brackets. The cleanest possible bracket debonding with no ceramic splinters to damage soft tissue.




709-231-130
Direct Bond Bracket Remover
 13 cm, 5"
 Heavy duty construction makes light work of removing direct bond brackets. Diamond honed, beveled contact points work equally well with metal and plastic brackets.



709-234-130
Posterior Band Remover
 13 cm, 5"
 The replaceable plastic bushing rests on occlusal surface. Inserted beak wedges under band. long style arm rests on 5mm bushing. May be autoclaved with bushing removed.



709-237-130
Angulated Direct Bond
Remover
 13 cm, 5"
 For difficult work in confined intraoral areas provide significantly less patient discomfort. Choose from two models-both provides easy access to the far corners of the mouth. May be autoclaved with replaceable bushing removed.



709-240-130
Mid-sized Bird Beak Plier
(light Wire)
 13 cm, 5"
 light wire forming plier is ideal for configuring loops, helices and hooks in wires up to .028". Elimination of cutter allows for larger jaw platform. Has groove in tip of flat beak.

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



709-243-130
Flat Nose Plier
 13 cm, 5"



709-246-130
Tweed Arch Bending Plier
 13 cm, 5"
 use this specialty plier to form offset and 90 bends on rectangular and square wires up to .022" x .028". Radiuses on tip edges do not mar or fray wire during torquing.



709-246-130
Loop Forming Plier
 13 cm, 5"
 Ideal for making consistent omega loops in round and rectangular wires. Each replaceable tip has two diameters .045" and .060"



709-249-130
Nitanium Expander adjusting plier
 13 cm, 5"
 specifically designed for making precision adjustments to the titanium palatal expander and molar rotator. The grooved beaks assure a strong, non-slip hold during adjustment.



709-252-130
Flush Cut Distal End Cutters pliers
 13 cm, 5"



709-255-130
Nance Closing Loop Plier
 13 cm, 5"



709-258-130
Nickel Titanium Plier
 13 cm, 5"



709-261-130
De La Rosa Plier
 13 cm, 5"

Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



709-264-130
Direct Bonding Bracket
Remover (wide) Plier T.c
 13 cm, 5"



709-267-130
Hook Crimping Plier
 13 cm, 5"




709-270-130
Tube/stop Crimpable
Plier
 13 cm, 5"



709-273-130
Standard Angle Distal
End Flush Cutter
 13 cm, 5"
 Cuts round and multi-strand wire flush to the buccal tube without fraying the wire. A cotton roll behind the buccal tube when cutting helps to catch the end of the wire.



709-276-130
Reverse Action Distal
End Cutter
 13 cm, 5"
 The cutting edges to guarantee absolutely flush trims. Cuts round and multi-strand wires.



709-279-130
Safety Hold Distal End Cutter
 13 cm, 5"
 This best seller works on hard round, rectangular and braided wires. Safety hold eliminates flying wire shards. Cuts and holds in one operation.




709-282-130
Miniature Safety Hold
Distal End Cutter
 13 cm, 5"
 Compact design gives greater range of motion in tight quarters - with the same safety hold feature of standard size model.

*Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik*

*Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis*



709-285-130
Pin and ligature
cutter T.C
13 cm, 5"



709-288-130
Micro Ligature Cutter
13 cm, 5"
Overall dimensions are 10% smaller than our standard Pin and Ligature cutter. That means easier access in small mouths, and fatigue-free comfort for orthodontists with smaller hands. Our smallest model. New long handle version and a long handle with a 15 angled head.



709-291-130
Miniafture Pin & Ligature Cutter
13 cm, 5"
Compact head design allows easy access under bracket flange. Cut up 0.014". For elastic modules also.



709-293-130
Micro Miniafture Pin & Ligature Cutter
13 cm, 5"
compact size tip means easier access offers better control. Ideal for young patients. Cut up to .014"



709-296-130
Hard Wire Cutter
13 cm, 5"
hardened cutting edges make short work of hard round and rectangular wires. Slightly larger head size for increased durability, standard handle length for a sure, non-slip grip.



709-299-130
Distal End Cutter T.C
13 cm, 5"
Safety hold is ideal for cutting wires


PLIERS FOR ORTHODONTICS & PROSTHETICS

*Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik*

*Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis*



**709-302-130
Hard Wire Cutter**
13 cm, 5"
Straight



Cut Hard wire up to
.022" x .028"
(.56mm x .71mm)



**709-305-130
Pin & Ligature Cutter**
13 cm, 5"
The best of both worlds—small cutting head for better mobility and view of operative site plus full sized handles for better heft and feel. Tapered tips easily reach difficult areas.



**709-308-130
Oblique Pin & Ligature Cutter**
13 cm, 5"
By setting the cutting head at a 15 angle to the handles, you get better, more visible access to the working area. More professionals prefer oblique angled pin and Ligature cutters for wires up to .020".



**709-311-130
Distal End Cuter
Safety Holder**
13 cm, 5"
(Noninserted)



**709-314-130
Standard Pin & Ligature Cutter**
13 cm, 5"
Clean, flush cuts on pins, soft wire up to .020". Standard size plier body offers maximum comfort for large hands.



**709-317-130
Angular Ligature Cutter**
13 cm, 5"
Slim 15 nose design allows easy access under bracket tie wings for fast, sure and easy cutting of ligatures. Tapered profile and tip slant make flush cuts a snap.



**709-320-130
Angular Ligature Cutter**
13 cm, 5"
Tapered cutting edges set at 45 angle eases access to hard to reach intraoral areas. Designed to prevent cheek impingement.

*Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik*

*Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis*



709-323-130
Hook Crimping Pliers T.c
13 cm, 5"
angled straight



709-326-130
Distal End Cutter (non-holding).
13 cm, 5"
A precision cutter designed for distal arch cutting. Cutting angle of 72 degree allows cutting of wire very close to distal end of buccal tube. Caution: maximum cutting capacity arch wire up to 0.56x0.71 mm rectangular, round wires up to 0.51 mm.



709-329-130
Arch Forming Plier.
13 cm, 5"
Angled beak allows easy placement from the lingual of the anterior segment of edgewise of edge wires arch wires and increases forquing capability when used in a two-plier technique



709-332-130
Direct Bonding Bracket Remover (wide).
13 cm, 5"
The tips wedge between edges of brackets base and tooth surface, lifting bracket off with virtually no stress.



709-335-130
Tweed Arch Bending Plier.
13 cm, 5"
Modified with extra short 6.5mm and thin blades.



709-338-130
Aderer max Plier
13 cm, 5"
Ut. 0.7mm.



709-341-130
Utility Plier Weingart.
13 cm, 5"



709-344-130
Band Seating Plier.
13 cm, 5"



PLIERS FOR ORTHODONTICS & PROSTHETICS

*Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik*

*Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis*




**709-347-130
Utility Plier
Weingart T.c**
13 cm, 5"



**709-350-130
Ligual Arch Placing
and Removing Plier.**
13 cm, 5"



**709-353-130
Band Silting Plier**
13 cm, 5"
with stop screw.
Used to cut lower
and upper bands.



**709-356-130
Palatal Arch Plier**
13 cm, 5"
For 0.76mm & 90mm.



**709-359-130
Diagonal Cutter-
Compound Action**
13 cm, 5"
The compound action
cutter bends and cuts
wires up to 1.30mm (.051").



**709-362-130
How Straight
Plier**
13 cm, 5"



**709-365-130
Lingual Arch Plier**
13 cm, 5"
used to form double
and triple back bends
in 0.76 or 0.91 wire
for insertion in lingual
sheaths.



**709-368-130
Adams Bending Plier**
13 cm, 5"
for wire up to 0.76mm.
Lenght of the plier 11.5 cm



**709-371-130
Tweed Loop Forming Plier.**
13 cm, 5"
For making loops in square
or rectangular wires up
to 0.56x0.71mm. The
replaceable beak has three
diameter sections:
1.14mm-1.52mm-1.90mm.



Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik

Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis



709-374-130
V Bend Stop Plier T.c
 13 cm, 5"



709-377-130
Reynolds Contouring Plier.
 13 cm, 5"




709-380-130
Single Diameter Loop Farming Plier.
 13 cm, 5"
 The replaceable cylindrical back has only one diameter section of 1.52mm.



709-383-130
Omega Loop Forming Plier.
 13 cm, 5"
 For making perfect omega loops in square or rectangular wires up to 0.56x0.71mm.



709-386-130
45 Offset Angled Tips Plier.
 13 cm, 5"



709-389-130
Aderer max Plier
 13 cm, 5"
 Ut. 0.7mm.



709-392-130
Arch Forming Plier
 13 cm, 5"
 non grooved. Modified with a non grooved contouring surface for the forming of precise arches.



709-395-130
De La Rosa Arch Forming Plier
 13 cm, 5"
 grooved. Designed to form and contour arch wires with or without torquing on round or rectangular wire. The beak has grooves of 0.41mm, 0.46mm, 0.56mm

*Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik*

*Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis*



709-398-130
Direct Bond Bracket
Remover T.c
13 cm, 5"
Straight



709-401-130
Posterior Band
Remover Plier
13 cm, 5"
With Pad long Type.



709-404-130
Direct Bonding Bracket
Remover Plier
13 cm, 5"
with pad. Narrow point.



709-407-130
Direct Bonding Bracket
Remover Plier
13 cm, 5"




709-410-130
Double Rounded Jaw
Wire Bending Plier.
13 cm, 5"
Wires up to 0.76mm.



709-413-130
Round And Concave Plier
13 cm, 5"
for loop forming and bending
of thin round and edgewise
wires up to 0.51mm.



709-416-130
Light Wire Bird
Beak Plier
13 cm, 5"
for bending round,
square or rectangular
wire up to 0.41mm.



709-419-130
Bird Beak Plier
13 cm, 5"
for bending round,
square or rectangular
wire up to 0.76mm.

*Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik*


*Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis*



**709-422-130
Standard Ligature
Cutter T.c**
13 cm, 5"



**709-425-130
Round Nose Wire
Bending Plier.**
13 cm, 5"
Wire up to 0.51mm.



**709-428-130
Flat Nose Plier**
13 cm, 5"
with smooth beaks
for bending wire up
to 0.76mm.



**709-431-130
Micro-Miniature
Light Wire Cutter**
13 cm, 5"
(with stop screw).
Capacity: soft ligature
wire up to 0.30mm.



**709-434-130
Ligature Cutter.**
13 cm, 5"
45° angle
ligature wires and pins.
Capacity: ligature wires
Tip to 0.38mm. At the tips.



**709-437-130
Tweed Arch Bending
Plier-Short.**
13 cm, 5"
Plier with short blades
10mm. Designed for
convenient handling of
square or rectangular
wire.



**709-440-130
Tweed Rectangular
Arch Forming Plier.**
13 cm, 5"
Plier with long
blades 14mm.



**709-443-130
Direct Bond Attachment
Remover (angled).**
13 cm, 5"
For removing all bonded
attachments. May be
used mesial-distal (horizontal)
or incisal-gingival (vertical) in
anterior areas.



**709-446-130
Distal End Cutter.**
13 cm, 5"
Safety hold is ideal for cutting
wires within the patient's mouth.
Plier shear-cuts hard wires
close to the buccal tube and
holds loose wire.(maximum)
cutting capacity: archwires up
to .022" x .028")

*Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik*

*Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis*



709-449-130
Mini Pin & Ligature Cutter T.c
13 cm, 5"



709-452-130
Small Distal End Cutter
13 cm, 5"
safety hold is ideal for cutting wires within the patient's mouth. Plier shear-cuts hard wires close to the buccal tube and holds loose wire. (Maximum cutting Capacity: archwires up to .022" x .028")



709-455-130
Lab Hard Wire Cutter/Chrome.
13 cm, 5"
The extra strong edges make it one of the most popular hard wire cutters. Designed to cut all sizes and types of archwires, it features a 15 cutting angle for easy intraoral access. (Cuts archwires up to .022" x .028")



709-458-130
Mini Pin & Ligature Cutter
13 cm, 5"
Fine tips to easily access hard to reach areas. Cuts wire ligatures, pins and elastics. Maximum cutting Capacity: soft wire up to .015". 7 cutting angle



709-461-130
Micro Mini Pin & Ligature Cutter
13 cm, 5"
Compact tip size allows easy access to tight inter-bracket areas. Super-sharp, diamond honed edges are effective for cutting all soft wires, pins and elastics up to .012". 7 cutting angle



709-464-130
Heavy Duty Wire Cutter
13 cm, 5"
Unique joint and extremely sturdy cutting edge make it ideal for cutting face bows and other large diameter wires up to .08". This cutter is not autoclavable.



709-467-130
Angled Direct Bond Bracket Removing Plier
13 cm, 5"
Designed to quickly and sagely remove direct bond brackets. The 60 angle allows access to anterior and posterior brackets

*Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik*

*Pinces pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis*



**709-470-130
Angled Direct
Bond Bracket
Remover T.c**
13 cm, 5"



**709-473-130
How Plier**
13 cm, 5"
Matching serrated pads assure a firm, positive grip. Rounded tips for patient comfort and safety. 3/32" diameter pads.



**709-476-130
Angle How Plier**
13 cm, 5"
Efficient utility plier with serrated beaks for a firm grip. 45° angled tips to access hard to reach areas.



**709-479-130
Weingart Plier**
13 cm, 5"
Serrated tips hold wires firmly at a convenient working angle. Tapered beaks fit easily between brackets and hard to reach areas. Our most popular utility plier.



**709-482-130
Light Wire Bending Plier**
13 cm, 5"
Designed for working on small wires and forming wire springs. Tapered beaks make it easy to bend small diameter loops. (Bends wires up to .020")



**709-485-130
Jarabak Plier**
13 cm, 5"
Excellent for precise bending and forming loops in light wires up to .020" 3 set of precision grooves assure a firm grip.



**709-488-130
Flat On Flat Plier**
13 cm, 5"
Stubby, square beaks; .045" at tip and parallel at 1mm opening. Easily bends and forms wire up to .050".



**709-491-130
Occulist Plier**
13 cm, 5"
Round and concave beaks bend round or rectangular wire with a firm grip that will not nick wire. Bends wires up to .020".



**709-494-130
Hollow chop arch
forming plier**
13 cm, 5"



**709-497-130
Tweed Loop Forming Plier**
13 cm, 5"
convex/concave beak makes precise omega loops. Ti has 3 forming sections: .045", .060". And .075". Etched tip surfaces prevent wire slippage.

*Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik*


*Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis*



709-500-130
3 - Prong Clasp
Adjustable Plier T.c
13 cm, 5"



709-503-130
Loop Forming & Closing Plier
13 cm, 5"
Precision loop forming with four step tip (3,4,5 & 6mm). Excellent for wires up to .022". Beveled edges prevent wire scoring & one serrated beak prevents wire slippage.



709-506-130
"V" Bend stop Plier
13 cm, 5"
Special design makes 1mm "V" Bends to shorten archwires or provide a positive stop. Excellent for placing stops in Nitanium wires. For round or rectangular wires up to .025".



709-509-130
Loop Tie Back Plier
13 cm, 5"
Four-stop plier automatically forms loops on wire up to .020" with one motion. Closing loops are easily made to desired height.



709-512-130
3-prong Clasp Adjusting Plier
13 cm, 5"
Precision tips for accurate wire and clasp adjusting and contouring. Also an excellent plier for bending of orthopedic appliances. Superior strength for wires up to .030".



709-515-130
Torquing Plier Set With Key
13 cm, 5"
Specially designed to place torquing bend into the small sections of archwires (such as a bend for a single tooth) without distorting adjacent sections of the wire. Two piece set consists of one torquing plier and one .018/.022 key.



709-518-130
Detailing Step Plier
13 cm, 5"
Double-sided offset beak permits bayonet bends. Creates right hand or left hand step up or step down at 1/4, 3/4, or 1 mm increments.



709-521-130
Hook Crimping Plier
13 cm, 5"
Designed for placement and securing of crimpable archwire hooks & power hooks. Angled beak allows for precise placement of hook on arch wire.

*Pliers for orthodontics and prosthetics
Zangen für Kieferorthopädie und Prothetik*

*Pincès pour l'orthodontie et la prothèse
Alicates para ortodoncia y prótesis*



**709-524-130
V Bend Stop
Plier T.c**
13 cm, 5"



**709-527-130
Pin Holding
Plier T.c**
13 cm, 5"



**709-530-130
Direct Bracket
Remover T.c**
13 cm, 5"



**709-533-130
Hook Crimping
Plier**
13 cm, 5"



**709-536-130
Angled Utility Arch Plier**
13 cm, 5"
Angled beak with pin simplifies the placement and manipulation of preformed utility archwires. Forms 3mm steps for accurate utility arch construction.



**709-539-130
Lingual Arch Forming Plier**
13 cm, 5"
Plier makes consistent, repeatable bends. Forms double-back and triple-back bends in .030" wire for use with lingual sheaths.



**709-542-130
E-Z Ligature Placer**
13 cm, 5"
this metal-tipped tool makes application of elastic ligatures fast and easy. Place closed plier tip into ligature, pull ligature from cane, stick or ring, then squeeze plier to stretch elastic open and apply ligature.

Orthodontics Instruments
Orthodontische Instrumente

Instrumentos Ortodonticos
Instrumentos Ortodontiques



709-545-001
Ortho Light Band Pusher Scaler



709-548-002
Ortho Light double ended ligature director



709-551-003
Ortho Light double ended ligature director scaler



709-554-004
Ortho Light double ended ligature scaler/director
tivate coil spring are remove modules with one instruments



709-557-005
Ortho band pusher scaler



709-560-006
Ortho Light Band Pusher Scaler



709-563-007
Ortho ligature director single ended



709-566-008
Ortho Light double ended ligature director scaler



709-569-009
Ortho Light double ended ligature scaler/director
tivate coil spring are remove modules with one instruments

Orthodontics Instruments
Orthodontische Instrumente

Instrumentos Ortodonticos
Instruments Orthodontiques



709-572-010
Ortho Light double ended ligature director



709-575-011
Ortholite Ligature Director Scaler



709-578-012
Double ended cement spatula



709-581-013
Ortho light ligature director double ended



709-584-014
Ortho light scaler single ended



709-587-015
Ortho light scaler double ended



709-590-016
Ortho light scaler double ended



709-593-017
Ortho light scaler double ended

Orthodontics Instruments
Orthodontische Instrumente

Instrumentos Ortodonticos
Instruments Orthodontiques



709-596-018
Twiel-on elastic module instrument



709-599-019
Twiel-on elastic module instrument



709-602-020
Schure Band Seater



709-605-021
Ortho light scaler



709-608-022
Ortho light scaler



709-611-023
Ortho light scaler heavy



709-614-024
Ortho instrument



709-617-025
Ortho instrument

Orthodontics Instruments
Orthodontische Instrumente

Instrumentos Ortodonticos
Instrumentos Ortodontiques



709-620-026 Band seater auto cleavable



709-623-027 Band seater auto cleavable



709-626-028 Ortho instrument



709-629-029 Heavy duty band pusher

709-631-030
Bone Bracket Gauge



709-634-031 Positive Positioner



709-637-032 Distal Bender



709-640-033 Polar Bear Ice Instruments

Impression Trays-Stainless Steel
Abdruckloeffel S.S.




Cubetas para impresiones de acero inox
Porte-empresintes en acier inoxydable

SUPER REGULAR PERFORATED RIM LOCK (Set of 10 Pcs)

709-643-001, Set of 6 Trays-No: Upper-2/L, 3/M, 4/S and Lower-2/L, 3/M, 4/S

709-646-002, Set of 8 Trays-No: Upper-1/XL, 2/L, 3/M, 4/S and Lower-1/XL, 2/L, 3/M, 4/S

709-649-003, Set of 10 Trays-No: Upper-1/XL, 2/L, 3/M, 4/S, 5/XS and Lower-1/XL, 2/L, 3/M, 4/S, 5/XS **DENTULAS**



REGULAR TRAYS PERFORATED

	A	B
#1-Upper	3.05	2.45
#2-Upper	2.87	2.33
#3-Upper	2.74	2.15
#4-Upper	2.87	2.16
#5-Upper	2.49	1.06
#1-Lower	2.90	2.44
#2-Lower	2.83	2.33
#3-Lower	2.88	2.14
#4-Lower	2.87	2.18
#5-Lower	2.75	2.00



All dimensions are in inches

SUPER REGULAR SOLID RIM LOCK (Set of 10 Pcs)

709-651-004, Set of 6 Trays-No: Upper-2/L, 3/M, 4/S and Lower-2/L, 3/M, 4/S

709-654-005, Set of 8 Trays-No: Upper-1/XL, 2/L, 3/M, 4/S and Lower-1/XL, 2/L, 3/M, 4/S

709-657-006, Set of 10 Trays-No: Upper-1/XL, 2/L, 3/M, 4/S, 5/XS and Lower-1/XL, 2/L, 3/M, 4/S, **DENTULAS**



REGULAR TRAYS SOLID


	A	B
#1-Upper	3.04	2.44
#2-Upper	2.84	2.27
#3-Upper	2.68	2.11
#4-Upper	2.54	2.19
#5-Upper	2.47	1.99
#1-Lower	2.94	2.41
#2-Lower	2.82	2.30
#3-Lower	2.44	2.19
#4-Lower	2.76	2.06
#5-Lower	2.80	1.09

All dimensions are in inches

IMPRESSION TRAYS-STAINLESS STEEL

*Impression Trays-Stainless Steel
Abdruckloffel S.S.*

*Cubetas para impresiones de acero inox
Porte-empresintes en acier inoxydable*




ORTHODONTIC (CHILDREN) PERFORATED RIM LOCK (Set of 8 Pcs)

709-660-007, Set of 8 Trays-No: Upper-15,16,17,18 and Lower-15,16,17,18

ORTHODONTIC TRAYS PERFORATED

	A	B
#15-Upper	2.49	2.00
#16-Upper	2.44	1.94
#17-Upper	2.41	1.71
#18-Upper	2.13	1.69
#15-Lower	2.78	2.08
#16-Lower	2.73	1.06
#17-Lower	2.52	1.92
#18-Lower	2.24	1.65

All dimensions are in inches




ORTHODONTIC (CHILDREN) SOLID RIM LOCK (Set of 8 Pcs)

709-663-008, Set of 8 Trays-No: Upper-15,16,17,18 and Lower-15,16,17,18

ORTHODONTIC TRAYS SOLID

	A	B
#15-Upper	2.48	2.03
#16-Upper	2.49	1.94
#17-Upper	2.38	1.68
#18-Upper	2.12	1.70
#15-Lower	2.77	2.08
#16-Lower	2.75	1.95
#17-Lower	2.45	1.95
#18-Lower	2.30	1.69

All dimensions are in inches





*Impression Trays-Stainless Steel
Abdruckloeffel S.S.*

*Cubetas para impresiones de acero inox
Porte-empreintes en acier inoxydable*

EDENTULOUS (FULL DENTURE) PERFORATED RIM LOCK (Set of 8 Pcs)

709-666-009, Set of 6 Trays-No: Upper-2/L, 3/M, 4/S and Lower-6/L, 7/M, 8/S
709-669-010, Set of 8 Trays-No: Upper-1/XL, 2/L, 3/M, 4/S and Lower-5/XL, 6/L, 7/M, 8/S

EDENTULAS




TRAYS PERFORATED		
	A	B
#1-Upper	3.03	2.68
#2-Upper	2.84	2.40
#3-Upper	2.62	2.50
#4-Upper	2.58	2.29
#5-Lower	3.27	2.44
#6-Lower	2.97	2.40
#7-Lower	2.85	2.50
#8-Lower	2.57	2.20

All dimensions are in inches

EDENTULOUS (FULL DENTURE) SOLID RIM LOCK (Set of 8 Pcs)

709-672-011, Set of 6 Trays-No: Upper, 2/L, 3/M, 4/S and Lower-6/L, 7/M, 8/S
709-675-012, Set of 8 Trays-No: Upper-1/XL, 2/L, 3/M, 4/S and Lower-5/XL, 6/L, 7/M, 8/S

EDENTULAS



TRAYS SOLID		
	A	B
#1-Upper	3.01	2.67
#2-Upper	2.88	2.28
#3-Upper	2.57	2.42
#4-Upper	2.53	2.33
#5-Lower	3.23	2.40
#6-Lower	3.00	2.35
#7-Lower	2.85	2.14
#8-Lower	2.62	2.22

All dimensions are in inches

IMPRESSION TRAYS FOR PARTIAL IMPRESSIONS S.S.



*Impression Trays For Partial Impressions S.S.
Abdruckloeffel fur partielle abdrucke S.S.*

*Cubetas para impresiones parciales S.S.
Porte-empresiones Pour empreintes Particelles S.S.*

PARTIAL PERFORATED RIM LOCK

- 709-678-013, Partial Upper Left and Lower Right
- 709-681-014, Partial Upper Right and Lower Left
- 709-684-015, Partial Anterior Lower
- 709-687-016, Partial Anterior Upper
- 709-690-017, Partial Turnable

- 709-693-018, Set of 2- #S. 30,31
- 709-696-019, Set of 2- #S. 32, 33
- 709-699-020, Set of 3- #S. 30, 31, 99
- 709-700-021, Set of 4- #S. 30, 31, 32, 33
- 709-703-022, Set of 5- #S. 30, 31, 32, 33,99





PARTIAL TRAYS PERFORATED		
	A	B
#30-Partial Upper Left and Lower Right	1.40	2.45
#31-Partial Upper Right and Lower Left	1.40	2.45
#32-Partial Interior Lower	2.57	1.50
#33-Partial Interior Upper	2.57	1.50
#99-Partial Turnable	2.47	1.12

All dimensions are in inches

PARTIAL SOLID RIM LOCK

- 709-706-023, Partial Upper Left and Lower Right
- 709-709-024, Partial Upper Right and Lower Left
- 709-712-025, Partial Anterior Lower
- 709-715-026, Partial Anterior Upper
- 709-718-027, Partial Turnable

- 709-721-028, Set of 2- #S. 30,31
- 709-724-029, Set of 2- #S. 32, 33
- 709-727-030, Set of 3- #S. 30, 31, 99
- 709-730-031, Set of 4- #S. 30, 31, 32, 33
- 709-733-032, Set of 5- #S. 30, 31, 32, 33,99



PARTIAL TRAYS SOLID		
	A	B
#30-Partial Upper Left and Lower Right	1.40	2.45
#31-Partial Upper Right and Lower Left	1.40	2.45
#32-Partial Interior Lower	2.57	1.50
#33-Partial Interior Upper	2.57	1.50
#99-Partial Turnable	2.47	1.12

All dimensions are in inches

IMPRESSION TRAYS-STAINLESS STEEL

Impression Trays-Stainless Steel
Abdruckloffel S.S.



Cubetas para impresiones de acero inox
Porte-empreintes en acier inoxydable

Water Cooled Solid Impression Trays for Hydrocalliod Impression Complete Set.
5 - Upper, 5 - Lower
709-736-033, Set of 10 Trays

Check Bite
709-739-034


Upper, 1/XL, 2/L, 3/M, 4/S, 5/XS

Lower, 1/XL, 2/L, 3/M, 4/S, 5/XS



REGULAR SOLID RIM LOCK


709-742-035, Set of 6 Trays No:
Upper - 2/L, 3/M, 4/S and Lower - 2/L, 3/M, 4/S
709-745-036, Set of 8 Trays No:
Upper - 1/XL, 2/L, 3/M, 4/S and
Lower - 1/XL, 2/L, 3/M, 4/S
709-748-037, Set of 10 Trays No:
Upper-1/XL, 2/L, 3/M, 4/S, 5/XS
and Lower - 1/XL, 2/L, 3/M, 4/S, 5/XS
709-751-038, Set of 8 Trays No:
Upper-2/L, 3/M, 4/S, 5/XS
and Lower - 2/L, 3/M, 4/S, 5/XS




REGULAR PERFORATED RIM LOCK

709-754-039, Set of 6 Trays No:
Upper - 2/L, 3/M, 4/S and Lower - 2/L, 3/M, 4/S
709-757-040, Set of 8 Trays No:
Upper - 1/XL, 2/L, 3/M, 4/S and
Lower - 1/XL, 2/L, 3/M, 4/S
709-760-041, Set of 10 Trays No:
Upper-1/XL, 2/L, 3/M, 4/S, 5/XS
and Lower - 1/XL, 2/L, 3/M, 4/S, 5/XS
709-763-042, Set of 8 Trays No:
Upper-2/L, 3/M, 4/S, 5/XS
and Lower - 2/L, 3/M, 4/S, 5/XS




709-766-001Complete Instrument Tray With Perforated Base,
Lid & Insert Frame For 16 Instruments**709-772-003**


Lid Dimensions 288mm x 187mm x 29mm

**709-778-005**


Professional Base: Dimension 284mm x 183mm x 17mm

709-769-002Complete Instrument Tray With Perforated Base,
Lid & Insert Frame For 23 Instruments**709-775-004**


Lid Dimensions 288mm x 187mm x 29mm

**709-781-006**


Professional Base: Dimension 284mm x 183mm x 17mm

709-784-007Complete Instruments Tray With Solid Base,
Lid & Insert Frame For 16 Instruments**709-790-009**


Lid Dimensions 288mm x 187mm x 29mm

**709-796-011**

Solid Base: Dimensions 284mm x 183mm x 17mm

709-787-008Complete Instruments Tray With Solid Base,
Lid & Insert Frame For 23 Instruments**709-793-010**


Lid Dimensions 288mm x 187mm x 29mm

**709-799-012**

Solid Base: Dimensions 284mm x 183mm x 17mm



709-802-015
6" x 8"
INSTRUMENT TRAY
With Lid




709-805-014
6" x 10"
INSTRUMENT BOX




709-808-017
INSTRUMENT TRAY
With Tray



709-811-018 6" x 8"
709-814-019 8" x 10"
709-817-020 10" x 12"
Without Lid




709-820-023
ENDODONTIC BOX
Also Available in Solid Pattern




709-823-024
8" x 12"
With Lid



NEEDLE CASE
709-830-026
230 x 125 mm



SCALER BOX
709-833-027
300 x 180 mm



709-836-028
Scallers Tray



SCALER TRAY
709-839-029
100 x 200 mm



INSTRUMENT TRAY WITHOUT COVER

	Size		Size
709-842-030	240 x 240 mm	709-875-041	230 x 125 mm
709-845-031	355 x 400 mm	709-878-042	200 x 150 mm
709-848-032	300 x 150 mm	709-881-043	250 x 200 mm
709-851-033	300 x 180 mm	709-884-044	300 x 250 mm
709-854-034	300 x 240 mm	709-887-045	350 x 300 mm
709-857-035	340 x 320 mm	709-889-046	400 x 350 mm
709-860-036	230 x 180 mm	709-892-047	450 x 300 mm
709-863-037	250 x 130 mm	709-895-048	300 x 200 mm
709-866-038	200 x 130 mm	709-898-049	350 x 250 mm
709-869-039	180 x 130 mm	709-901-050	250 x 125 mm
709-872-040	150 x 100 mm	709-904-051	200 x 100 mm
		709-907-052	400 x 200 mm



KIDNEY TRAY WITH COVER
(English Pattern)

	Size
709-910-053	150 mm
709-913-054	200 mm
709-916-055	240 mm
709-919-056	300 mm



KIDNEY TRAY (German Pattern)

	Size
709-922-057	150 mm
709-925-058	200 mm
709-928-059	240 mm
709-931-060	300 mm



STERILIZING DRUM

	Size
709-934-061	150 x 150 mm
709-937-062	175 x 175 mm
709-940-063	190 x 190 mm
709-943-064	200 x 140 mm
709-946-065	200 x 200 mm
709-949-066	240 x 120 mm
709-952-067	225 x 225 mm
709-955-068	240 x 240 mm
709-958-069	290 x 160 mm
709-961-070	280 x 240 mm
709-964-071	230 x 150 mm
709-967-072	300 x 250 mm
709-970-073	340 x 250 mm
709-973-074	350 x 250 mm
709-976-075	390 x 200 mm



FORCEPS JAR

	Size
709-979-076	55 x 100 mm
709-982-077	55 x 140 mm
709-985-078	55 x 170 mm
709-988-079	55 x 180 mm
709-991-080	80 x 100 mm



FORCEPS JAR

	Size
709-994-081	55 x 100 mm
709-997-082	55 x 140 mm
709-1000-083	55 x 170 mm
709-1003-084	55 x 180 mm
709-1006-085	80 x 100 mm



709-1009-086
COTTON HOIDER



ROUND BOWLS
709-1012-087
80 x 35mm
125 x 60mm



709-1015-088
SPRIT LAMP



709-1018-089
SPRIT LAMP



BOWL

	Size
709-1021-090	80 mm
709-1024-091	100 mm
709-1027-092	125 mm
709-1030-093	150 mm
709-1033-094	200 mm
709-1036-095	250 mm
709-1039-096	300 mm
709-1042-097	350 mm
709-1045-098	400 mm



WASH BASIN

	Size
709-1048-099	80 mm
709-1051-100	100 mm
709-1054-101	125 mm
709-1057-102	150 mm
709-1060-103	200 mm
709-1063-104	250 mm
709-1066-105	300 mm
709-1069-106	350 mm
709-1070-107	400 mm



BOWL WITH LIP


	Size
709-1075-108	80 mm
709-1078-109	100 mm
709-1081-110	125 mm
709-1084-111	150 mm
709-1087-112	200 mm
709-1090-113	250 mm
709-1093-114	300 mm
709-1096-115	350 mm
709-1099-116	400 mm




709-1100-001
Dental X-Ray Hanger
Single Piece Hanger



709-1101-002
Dental X-Ray Hanger
6"
use in six film sizes
stainless steel



709-1102-003
Film Holder For
10 Films
709-1103-004
14 Films



709-1104-005
Napkin Holder
with plastic band



709-1105-006
Napkin Holder
with Metal Chain